

The iHelp Mini can bring you and your loved ones the peace of mind that comes from knowing you can reach emergency personnel quickly and easily if and when you need it.

It's built-in 2-way microphone allows you to both talk and listen to trained emergency monitoring center operators. Its GPS tracking system ensures you can always be located. And it's small, lightweight design means you can easily take it wherever you go.

# ihelp<sup>TM</sup> mini

*A reliable and easy-to-use device for getting medical help when you need it.*



Scan the QR code for product information and support.

For more information on the iHelp Mini, call your dealer, or go to:  
[www.iHelpAlarm.com](http://www.iHelpAlarm.com)

# Features & Functions

Meet the **iHelp Mini**. An advanced mobile, medical alert system that can help you remain active, independent and safe whether at home or on the go.

*Reliable. Adaptable. Intelligent.*



## Crystal Clear Audio

The 2-way speaker allows for both talking and listening with emergency personnel.

## Two-way Communication

Operates on various 4G bands, and will fall back to 3G/2G when there is no 4G coverage.



## Geo Alert

Allows for multiple Geo fences and sends alerts as the device enters or leaves a particular area.

## Power On/Off Alert

Voice alerts will notify you if the device is turned off and when it is turned on.

## What's Included

- iHelp Mini
- Magnetic Charging Cord
- AC Adapter
- Breakaway Lanyard



## Superior Construction

- Comfortable to wear
- Lightweight - 1.41 ounces
- Small size - 2.4" x 1.732" x 0.63"



## SOS Alarm

An SOS alarm with location details will be sent to a nearby monitoring station when you press the SOS button or automatic fall detection is activated.

## Fall Alarm

In the event of a fall or collision the device will automatically send an SOS message to a nearby emergency center.

## Low Battery Alarm

When the battery charge falls below 20% the device will emit a verbal "low battery" message.

## Preset Voice Prompts

You can set voice prompts to remind yourself of important times (e.g., medication reminder).

## Options for Wearing your iHelp Mini

Can be worn around the neck (breakaway lanyard included) or an optional belt clip or wristband can be purchased separately.

## Optional Accessories

- Charging Base
- Belt Clip
- Wristband